MGN0932

Valid from / en cours de validité depuis le / en espera de validación desde el / in corso di validità dal / gültig seit dem

20/01/2019

Product Sheet (EN) 2

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
Eigenschaften (DE) 10

PRODUCT SHEET

valid from 20/01/2019

IDENTIFICATION OF THE MEDICAL DEVICE	Type	Rechargeable battery	
	Commercial designation	Batterie médicale BCI Mini-Torr Plus 6V 2.5Ah FCI	
	Reference	MGN0932	
	EAN	3660766621047	
	Brand	NX	
	Compatible / Original battery	Compatible	
	Packaging	Unitary	

RECOMMENDED USAGE

Follow the instructions and recommendations specific to each model, using the technical instructions and document resources from the devices in which the battery is used

Brands	Equipment	Models	PN
Smiths Medical	Monitor	Biochem Mini-Torr Plus	58522B2

Identification	GENERAL TECHNICAL CHARACTERISTICS	Chemistry	NiCd
		Type	C
		IEC designation	5 KRMU 26/50
		Rated voltage	6 V
		Nominal capacity	2.5 Ah
		Internal resistance Ω	< 40m Ω

The voltage and the actual capacity in use can be affected by several factors, especially the temperature, the discharge current, the pack's history (ex:use, storage), etc

ELECTRICAL CHARACTERISTICS	CHARGE	Maximum charging current	8.5 V
		Standard charging current (15h)	250 mA
		Fast charging current (2,5h)	N/A
	DISCHARGE	Range of operating voltage	5.5V at 8.5V
		Min tension in discharge	5V
		Max discharge current	3A
		Lifespan 80% DOD (0,5 C)	> 500 Cycles
	MAINTENANCE	Frequency of maintenance charges at 20°C	6 Months
	CONTROL ELECTRONICS	Electrical protection	Yes
		Low voltage power cut	No
High voltage power cut		No	
Max power cut voltage		Yes	

These devices not only designed to protect the pack in case of an equipment failure. They must not be used to control the discharge. The protection circuits have a response time of a few milliseconds.

MECHANICAL CHARACTERISTICS	Dimensions (+/- 2mm)	Length	126 mm
		Width	50 mm
		Depth	26 mm
	Weight (+/- 5g)	350 Gr.	
	Mechanical protection	Sleeve	
	Wire length (+/- 10mm)	150mm	
	Terminal	FCI 78211	

CONDITIONS OF USE, STORAGE, AND TRANSPORT	CONDITIONS OF USE	Charging temperature	5 at 55°C
		Discharge temperature	5 at 55°C
	CONDITIONS OF STORAGE	Storage temperature	5 at 25°C
		Level of humidity	65.00 %
		Max storage time	2 Years
	TRANSPORT	UN code	2800
		ADR/RID classification	Classe 8
		IMDG classification	Classe 8
		IATA classification	Classe 8

INSTRUCTIONS	COMMISSIONING	<ul style="list-style-type: none"> • Check the batteries and the connectors: wires not damaged, battery not swollen, burnt smell, oxidation of the connectors, leak... • Respect the polarity • Do a full charge with the adequate charger before the first use
	CHARGE	<ul style="list-style-type: none"> • Use an adequate charger • The battery is warmer during the charge: during the first charge, check that the battery's temperature stays in the temperature operating ranges. • In case of an abnormal heating, stop the charge by unplugging the charger within the realms of possibility, remove the battery from the equipment de l'équipement, have the equipment, the battery and the charger checked by a technician.
	CASE OF NON-WATERPROOF BATTERIES	<ul style="list-style-type: none"> • It is normal to observe a release of gas during the charge and use. Do not smoke. Place in suitable premises. • Open batteries need regular maintenance carried out by a qualified technician.
	CASE OF LITHIUM ION BATTERIES	There is a fire hazard with lithium ion batteries in the following cases: overload, short circuit, charge and use outside the voltage and temperature ranges.
	WARNINGS	<ul style="list-style-type: none"> • Read the instructions of your device. • Only use in compatible devices. • Respect the load and storage conditions. • Do not use if the battery is damaged, do not burn, do not pierce, do not dismantle or modify. The protection circuits protect the battery and the equipment: do not deactivate them.

Marking	CE
----------------	-----------

EXPLANATION OF SYMBOLS	LABELLING	
		Catalogue reference
		Lot number
		manufacturer's address
		Conform to the directive CE 93/42 DM class I
		To recycle in a suitable salvage and recycling structure
	Read the product sheet and the instruction manual	

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