



Rechargeable medical battery

Batterie médicale rechargeable

Batería médica recargable

MAIN INFORMATION | INFORMATIONS GÉNÉRALES | INFORMACIÓN GENERAL



Brand | Marque | Marca **CS**

Type of battery | Type de batterie | Tipo de batería

Compatible

Application | Application | Aplicación

Monitor / Moniteur / Monitor

Technology | Technologie | Tecnología

Nickel-métal hydrure

Voltage | Tension | Tensión

8,4V

Capacity | Capacité | Capacidad / Voltaje

8000mAh

Energy | Énergie | Energía

67200Wh

*Generic discharge value for the battery chemistry (not model-specific) | *Valeur générique de décharge propre à la technologie (non spécifique au modèle) | *Valor genérico de descarga asociado a la tecnología (no específico del modelo)

← DIMENSIONS | DIMENSIONS | DIMENSIONES →

Length | Longueur | Longitud **117mm** Width | Largeur | Anchura **69mm** Height | Hauteur | Altura **34mm**



Weight Battery
Poids batterie | Peso
0,85kg



Ean code | Code Ean
Codigo Ean
3660766629913

COMPATIBLE PRODUCTS* / PRODUITS COMPATIBLES* / PRODUCTOS COMPATIBLES*

PART NUMBER | PART NUMBER | PART NUMBER

2023227-001 **AMED2250** **M2828NO**
2023852-029 **M2828** **N1082**

FIT MODEL | MARQUES MODÈLES | MARCAS Y MODELAS

GE Healthcare / Datex Ohmeda / AMED2250

GE Healthcare / Datex Ohmeda / Physiologic Monitoring System

GE Healthcare / Datex Ohmeda / Dash 2500

Hellige / Marquette / Dash 2500

GE Healthcare / Datex Ohmeda / Dash 2500 Patient Monitor

Hellige / Hellige / Marquette Healthcare / Datex O / Hellige / Marquette MD

2500

GE Healthcare / Datex Ohmeda / Moniteur Dash 2500

GE Healthcare / Datex Ohmeda / Monitor Dash 2500

OUTLINE SAFETY WARNING: USE ONLY WITH IN THE ALLOWED PARAMETERS. Do not short circuit or over-load the battery. Charge only using an approved charger designed specifically to charge this battery. Do not heat above maximum temperatures indicated. Never crush, mutilate, puncture or abuse the battery. Do not dismantle the pack or disable any of the protective devices or circuits. DO NOT USE THE BATTERY IF YOU SUSPECT IT MAY BE FAULTY OR DAMAGED. © Copyright Enix Energies 2005. NB: This document and the product design are the intellectual property of Enix Energies. No document or design may be copied or used for commercial purposes without written permission of Enix Energies. Users must satisfy themselves, by means of testing etc., that products are suitable for their application. Data in this document is for general guidance only; consult cell manufacturers data for definitive information. Information is given free of charge and in good faith, but no responsibility can be accepted for any errors or omission or costs or losses or liabilities arising from the use of this information. All business is conducted to Enix Energies terms and conditions only.

