

Rechargeable medical battery
Batterie médicale rechargeable
Batería médica recargable

MAIN INFORMATION | INFORMATIONS GÉNÉRALES | INFORMACIÓN GENERALES



Brand Marque Marca	Arjo
Type of battery Type de batterie Tipo de batería	Original / Origine / Originale
Application Application Aplicación	Hoist / Lève-malade / Elevador de pacientes
Technology Technologie Tecnología	Nickel-métal hydrure
Voltage Tension Tensión	24V
Capacity Capacité Capacidad / Voltaje	2,5Ah
Energy Énergie Energía	60Wh

*Generic discharge value for the battery chemistry (not model-specific) | *Valeur générique de décharge propre à la technologie (non spécifique au modèle) | *Valor genérico de descarga asociado a la tecnología (no específico del modelo)

← DIMENSIONS DIMENSIONS DIMENSIONES →				
Length Longueur Longitud	Width Largeur Anchura	Height Hauteur Altura		
315mm	130mm	30mm	Weight Battery Poids batterie Peso 0,5kg	Ean code Code Ean Codigo Ean 3660766476463

COMPATIBLE PRODUCTS* / PRODUITS COMPATIBLES* / PRODUCTOS COMPATIBLES*

PART NUMBER | PART NUMBER |PART NUMBER

NEA0100
NEA0100-083

NEA0100083

FIT MODEL | MARQUES MODÈLES | MARCAS Y MODELAS

Arjo / Carendo	Arjo / BOB20 XX-01
Arjo / Carevo	Arjo / Carendo
Arjo / Carino	Arjo / Carevo
Arjo / Chaise douche BOB20 XX-01	Arjo / Carino
Arjo / Lève malade Mawi Twin	Arjo / Maxi twin
Arjo / Maxi twin	Arjo / Maxi twin compact
Arjo / Maxi twin compact	Arjo / Sara Combilizer
Arjo / Sara Combilizer	

OUTLINE SAFETY WARNING: USE ONLY WITH IN THE ALLOWED PARAMETERS. Do not short circuit or over-load the battery. Charge only using an approved charger designed specifically to charge this battery. Do not heat above maximum temperatures indicated. Never crush, mutilate, puncture or abuse the battery. Do not dismantle the pack or disable any of the protective devices or circuits. DO NOT USE THE BATTERY IF YOU SUSPECT IT MAY BE FAULTY OR DAMAGED. No document or design may be copied or used for commercial purposes without written permission of Enix Energies. Users must satisfy themselves, by means of testing etc., that products are suitable for their application. Data in this document is for general guidance only; consult cell manufacturers data for definitive information. Information is given free of charge and in good faith, but no responsibility can be accepted for any errors or omission or costs or losses or liabilities arising from the use of this information. All business is conducted to Enix Energies terms and conditions only.