

# Rechargeable medical battery

## Batterie médicale rechargeable

## Batería médica recargable

MAIN INFORMATION | INFORMATIONS GÉNÉRALES | INFORMACIÓN GENERAL



Brand   Marque   Marca	<b>Welch Allyn</b>
Type of battery   Type de batterie   Tipo de batería	<b>Original / Origine / Originale</b>
Application   Application   Aplicación	<b>Monitor / Moniteur / Monitor</b>
Technology   Technologie   Tecnología	<b>Lithium-ion</b>
Voltage   Tension   Tensión	<b>8,4V</b>
Capacity   Capacité   Capacidad / Voltaje	<b>1,9Ah</b>
Energy   Énergie   Energía	<b>15,96Wh</b>

\*Generic discharge value for the battery chemistry (not model-specific) | \*Valeur générique de décharge propre à la technologie (non spécifique au modèle) | \*Valor genérico de descarga asociado a la tecnología (no específico del modelo)

<p>← DIMENSIONS   DIMENSIONS   DIMENSIONES →</p> <p>Length   Longueur   Longitud    Width   Largeur   Anchura    Height   Hauteur   Altura</p> <p><b>128mm</b>                                <b>72mm</b>                                <b>35mm</b></p>	 <p>Weight Battery Poids batterie   Peso</p> <p><b>0,145kg</b></p>	 <p>Ean code   Code Ean Codigo Ean</p> <p><b>3660766476852</b></p>
--	---	---

COMPATIBLE PRODUCTS\* / PRODUITS COMPATIBLES\* / PRODUCTOS COMPATIBLES\*

PART NUMBER | PART NUMBER | PART NUMBER

**8085700**                                **008-0857-00**  
**10023002**                                **010-0230-02**

FIT MODEL | MARQUES MODÈLES | MARCAS Y MODELAS

- Welch Allyn / Micropaq**
- Welch Allyn / Propaq LT**
- Welch Allyn / Protocol Propaq 802LTRN Monitor**
- Welch Allyn / Micropaq**
- Welch Allyn / Propaq LT**
- Welch Allyn / Protocol Propaq 802LTRN Monitor**
- Welch Allyn /**

OUTLINE SAFETY WARNING: USE ONLY WITH IN THE ALLOWED PARAMETERS. Do not short circuit or over-load the battery. Charge only using an approved charger designed specifically to charge this battery. Do not heat above maximum temperatures indicated. Never crush, mutilate, puncture or abuse the battery. Do not dismantle the pack or disable any of the protective devices or circuits. DO NOT USE THE BATTERY IF YOU SUSPECT IT MAY BE FAULTY OR DAMAGED. No document or design may be copied or used for commercial purposes without written permission of Enix Energies. Users must satisfy themselves, by means of testing etc., that products are suitable for their application. Data in this document is for general guidance only; consult cell manufacturers data for definitive information. Information is given free of charge and in good faith, but no responsibility can be accepted for any errors or omission or costs or losses or liabilities arising from the use of this information. All business is conducted to Enix Energies terms and conditions only.